

## AUTOMATED ELECTROPHORESIS

# Experion<sup>™</sup> RNA HighSens and RNA StdSens Analysis Kits

- For automated electrophoresis, detection, and basic analysis with the Experion system
- Analysis of 1–12 samples in as little as 10–30 minutes
- Accurate RNA quantitation
- Single-step RNA purity determination
- Innovative Caliper Life Sciences'
   LabChip technology



## Single-Step RNA Analysis in 30 Minutes

#### Introduction

Experion RNA HighSens and RNA StdSens analysis kits combine state-of-the-art chip design with high-quality reagents to perform reproducible, quantitative, and accurate RNA analysis in a matter of minutes with the Experion system. Streamlined chip preparation methods and minimal sample requirements result in rapid experiments with minimal hands-on time.

# The Experion RNA HighSens and RNA StdSens Advantage

Innovative chip design, chemistry, software, and hardware make RNA separation and quantitation fast and easy.


- Single-step determination of RNA purity
- RNA ladder included in each kit
- Analysis of 1–12 samples in as little as 10–30 minutes
- Accurate quantitation of nanogram amounts of RNA

- Sensitivity down to 100 pg/µl with the RNA HighSens kit, dynamic range of 100–5,000 pg/µl
- Sensitivity down to 5 ng/µl with the RNA StdSens kit, dynamic range of 5–500 ng/µl

# **Experion RNA HighSens and RNA StdSens Chips**

Experion RNA chips offer an optimized microchannel architecture within an ergonomic caddy for exceptional RNA separation, chip handling, and results. Each RNA HighSens chip holds up to 11 samples, and each RNA StdSens chip holds up to 12. The chips are clearly labeled with well number, priming code, lot number, and expiration date for easy, error-free runs.


	Shelf Life	Storage Conditions
Experion RNA chips	9 months	Ambient
Experion RNA ladder	9 months	–20°C
Experion RNA reagents, supplies	9 months	4°C








## **Experience Meets Innovation**





#### Comparison of RNA separation on automated systems.

Upper panel, Experion RNA ladder separated on an RNA StdSens chip, displayed using Experion software; lower panel, a competitor's RNA ladder separated and displayed using a competitor's automated electrophoresis system. Run time, ~60 sec. The Experion RNA ladder provides more uniform peak heights (fluorescence intensities), resulting in improved quantitation.

## **Experion RNA HighSens and RNA StdSens Reagents**

Accurate assessment of RNA purity is essential to the success of quantitative RT-PCR and microarray experiments. Experion RNA analysis kits offer single-step RNA purity determination for confidence without delay. Experion RNA analysis kits include the following reagents and supplies:

- High-quality gel matrix for improved separation compared to agarose RNA gels, with significantly lower sample requirements
- Fluorescent dye for high-sensitivity, accurate detection
- Exclusive Experion RNA ladder for accurate sizing and quantitation
- Optimized loading buffer for reproducible results
- Sensitivity enhancer in RNA HighSens kit for accurate analysis at picogram levels
- Spin filters for convenient gel preparation

### **Ordering Information**

Catalog #	Description	Catalog #	Description
700-7103	Experion RNA StdSens Analysis Kit for 10 Chips, includes 10 RNA StdSens chips, 1,250 µl RNA gel, 20 µl RNA StdSens stain, 20 µl RNA ladder, 900 µl RNA StdSens loading buffer, 2 spin filters	700-7156	Experion RNA HighSens Reagents and Schips, includes 1,250 µl RNA gel, 20 µl stain, 20 µl RNA ladder, 900 µl RNA High buffer, 100 µl RNA sensitivity enhancer,
700-7104	Experion RNA StdSens Analysis Kit for 25 Chips, includes 25 RNA StdSens chips, 2 x 1,250 µl RNA gel, 2 x 20 µl RNA StdSens stain, 2 x 20 µl RNA ladder, 2 x 900 µl RNA StdSens loading buffer, 4 spin filters	700-7001	Experion System, 100/120 V, for RNA a includes electrophoresis station, priming station, software, USB2 cable, test chip (analysis kits sold separately)
700-7105	Experion RNA HighSens Analysis Kit for 10 Chips,	700-7252	Experion Electrode Cleaner, 250 ml
	includes 10 RNA HighSens chips, 1,250 µl RNA gel, 20 µl RNA HighSens stain, 20 µl RNA ladder, 900 µl RNA HighSens loading buffer, 100 µl RNA	700-7253  LabC	Experion DEPC-Treated Water, 100 ml Chip and the LabChip logo are trademarks
	sensitivity enhancer, 2 spin filters	" SpCUIP	Sciences, Inc. Bio-Rad Laboratories, Inc.
700-7106	Experion RNA HighSens Analysis Kit for 25 Chips, includes 25 RNA HighSens chips, 2 x 1,250 µl RNA gel,		aliper Life Sciences, Inc. to sell products logy for research use only.
	2 x 20 μl RNA HighSens stain, 20 μl RNA ladder, 2 x 900 μl RNA HighSens loading buffer, 2 x 100 μl RNA sensitivity enhancer, 4 spin filters	5,658,751, 5,43 and related forei	are licensed under U.S. Patent Nos. 5,86 6,134 and 5,582,977, and pending paten gn patents, for internal research and deve
700-7153	Experion RNA StdSens Chips, 10	, .	, quantitating and sizing macromolecules,
700-7154	Experion RNA StdSens Reagents and Supplies, for 10 chips, includes 1,250 µl RNA gel, 20 µl RNA StdSens stain, 20 µl RNA ladder, 900 µl RNA StdSens loading buffer, 2 spin filters	expressly excluding diagnostic, or ar	s, where internal research and developme des the use of this product for providing m ny other testing, analysis or screening serv I information or clinical analysis, in any eve
700-7155	Experion RNA HighSens Chips, 10	for compensatio	n by an unrelated party.

700-7156	Experion RNA HighSens Reagents and Supplies, for 10 chips, includes 1,250 µl RNA gel, 20 µl RNA HighSens stain, 20 µl RNA ladder, 900 µl RNA HighSens loading
	buffer, 100 µl RNA sensitivity enhancer, 2 spin filters
700-7001	Experion System, 100/120 V, for RNA analysis,
	includes electrophoresis station, priming station, vortex
	station, software, USB2 cable, test chip, instructions
	(analysis kits sold separately)
700-7252	Experion Electrode Cleaner, 250 ml

and the LabChip logo are trademarks of Caliper nces, Inc. Bio-Rad Laboratories, Inc. is licensed er Life Sciences, Inc. to sell products using the

licensed under U.S. Patent Nos. 5,863,753, 34 and 5,582,977, and pending patent applications, patents, for internal research and development use uantitating and sizing macromolecules, in combination here internal research and development use the use of this product for providing medical, ther testing, analysis or screening services, or ormation or clinical analysis, in any event in return y an unrelated party.

Practice of the patented polymerase chain reaction (PCR) process requires a license.



Bio-Rad Laboratories, Inc.

Life Science Group

Web site www.bio-rad.com USA (800) 4BIORAD Australia 02 9914 2800 Austria (01)-877 89 01 Belgium 09-385 55 11 Brazil 55 21 2527 3454 Canada (905) 712-2771 China (86 21) 6426 0808 Czech Republic + 420 2 41 43 05 32 Denmark 44 52 10 00 Finland 09 804 22 00
France 01 47 95 69 65 Germany 089 318 84-0 Greece 30 210 777 4396 Hong Kong (852) 2789 3300 Hungary 36 1 455 8800
India (91-124)-2398112/3/4, 5018111, 6450092/93 Israel 03 951 4127 Italy 39 02 216091 Japan 03-5811-6270 Korea 82-2-3473-4460
Latin America 305-894-5950 Mexico 55-52-00-05-20 The Netherlands 0318-540666 New Zealand 64 9 415 2280 Norway 23 38 41 30
Poland + 48 22 331 99 99 Portugal 351-21-472-7700 Russia 7 095 721 1404 Singapore 65-64153188 South Africa 00 27 11 4428508 Spain 34 91 590 52 00 Sweden 08 555 12700 Switzerland 061 717 95 55 Taiwan (886 2) 2578 7189/2578 7241 United Kingdom 020 8328 2000

Bulletin 3170 US/EG Rev B 04-0850 1204 Sig 1204